

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11108 - IMMUNIZING PEPTIDE

SKU: EBP11108

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_525128.2; NP_001035792.1 ID

 $\label{eq:lower_lower} \textbf{Immunogen} \ \, \begin{array}{l} \text{Peptide with sequence CRKKCRSQEYRIGR, from the internal region of the protein sequence according to NP_525128.2; NP_001035792.1.} \end{array}$

The peptide was used in the production of EB11108 - Goat Anti-DEFB104A / DEFB104B (aa37-50) Antibody. It is Product

Comments currently untested in blocking applications.

Peptide

CRKKCRSQEYRIGR Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade