



# EBP11109 - IMMUNIZING PEPTIDE

**SKU:** EBP11109



---

## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Shipped at ambient temperature, store at -20°C.
<b>Accession ID</b>	NP_689464.1; NP_001035794.1
<b>Immunogen</b>	Peptide with sequence CQKSLKCCRTIQP, from the internal region (near C Terminus) of the protein sequence according to NP_689464.1; NP_001035794.1.
<b>Product Comments</b>	The peptide was used in the production of EBP11109 - Goat Anti-DEFB106A / DEFB106B (aa47-59) Antibody. It is currently untested in blocking applications.
<b>Peptide Sequence</b>	CQKSLKCCRTIQP
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a>

## DOCUMENTS

- [Data Sheet](#)