

EBP11109 - IMMUNIZING PEPTIDE

SKU: EBP11109



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_689464.1; NP_001035794.1
Immunogen	Peptide with sequence CQKSLKCCRTIQP, from the internal region (near C Terminus) of the protein sequence according to NP_689464.1; NP_001035794.1.
Product Comments	The peptide was used in the production of EBP11109 - Goat Anti-DEFB106A / DEFB106B (aa47-59) Antibody. It is currently untested in blocking applications.
Peptide Sequence	CQKSLKCCRTIQP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)