

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11112 - IMMUNIZING PEPTIDE

SKU: EBP11112

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_695021.2; NP_775689.1; NP_697018.1 ID

 $\label{eq:lower_lower} \textbf{Immunogen} \ \ \ \ \ \ \ \ Peptide with sequence C-EENTDWSYEKQWPR, from the C Terminus of the protein sequence according to NP_695021.2; NP_775689.1; NP_697018.1.$

The peptide was used in the production of EB11112 - Goat Anti-Defensin, beta 119 Antibody. It is currently Product

Comments untested in blocking applications.

Peptide

C-EENTDWSYEKQWPR Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade