

EBP11123 - IMMUNIZING PEPTIDE

SKU: EBP11123



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_008872.1
Immunogen	Peptide with sequence C-QAHYKSAHLDHRHP, from the C or N Terminus of the protein sequence according to NP_008872.1.
Product Comments	The peptide was used in the production of EBP11123 - Goat Anti-SOX10 (aa204-217) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-QAHYKSAHLDHRHP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)