

EBP11136 - IMMUNIZING PEPTIDE

SKU: EBP11136



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_057419.4; NP_001182155.1
Immunogen	Peptide with sequence C-EDQESIQPLLRHYR, from the internal region of the protein sequence according to NP_057419.4; NP_001182155.1.
Product Comments	The peptide was used in the production of EBP11136 - Goat Anti-PRODH (aa112-134) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-EDQESIQPLLRHYR
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)