

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP11147 - IMMUNIZING PEPTIDE**

**SKU:** EBP11147

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

**Accession** 

ID

NP\_006763.2

Peptide with sequence C-ESAHIEREEYKLKEY, from the internal region (near C Terminus) of the protein sequence

according to NP\_006763.2.

The peptide was used in the production of EB11147 - Goat Anti-SYNGAP1 (aa1169-1183) Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

Sequence

C-ESAHIEREEYKLKEY

Shipping

Ambient Instructions

**Product** 

Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite\_plus\_medium.png$