

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11152 - IMMUNIZING PEPTIDE

SKU: EBP11152

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_149116.2 ID

 $\begin{tabular}{ll} \textbf{Immunogen} & Peptide with sequence ELSRCSRNKRENTS, from the internal region of the protein sequence according to NP_149116.2. \end{tabular}$

The peptide was used in the production of EB11152 - Goat Anti-SLC22A16 (aa97-110) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

ELSRCSRNKRENTS Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade