

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11159 - IMMUNIZING PEPTIDE

SKU: EBP11159

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 001185752.1; NP 031497.2; NP 001185753.1; NP 001185754.1; NP 001185755.1 ID

Peptide with sequence DAEFGHDSGFEVRHQK-C, from the N Terminus of the protein sequence according to Immunogen

NP 001185752.1; NP 031497.2; NP 001185753.1; NP 001185754.1; NP 001185755.1.

The peptide was used in the production of EB11159 - Goat Anti-Amyloid beta peptide (mouse) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

DAEFGHDSGFEVRHQK-C Sequence

Shipping

Grade

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png