

EBP11159 - IMMUNIZING PEPTIDE

SKU: EBP11159



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001185752.1; NP_031497.2; NP_001185753.1; NP_001185754.1; NP_001185755.1
Immunogen	Peptide with sequence DAEFGHDSGFVHRHQK-C, from the N Terminus of the protein sequence according to NP_001185752.1; NP_031497.2; NP_001185753.1; NP_001185754.1; NP_001185755.1.
Product Comments	The peptide was used in the production of EBP11159 - Goat Anti-Amyloid beta peptide (mouse) Antibody. It is currently untested in blocking applications.
Peptide Sequence	DAEFGHDSGFVHRHQK-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)