

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11173 - IMMUNIZING PEPTIDE

SKU: EBP11173

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

ID

NP_001129531.1

 $\begin{tabular}{l|l} \textbf{Immunogen} & Peptide with sequence C-QDRKLDENANVQLSD, from the internal region of the protein sequence according to NP_001129531.1. \end{tabular}$

Product The peptide was used in the production of EB11173 - Goat Anti-cytochrome P450 1A1 (mouse) Antibody. It is

Comments currently untested in blocking applications.

Peptide

Sequence

C-QDRKLDENANVQLSD

Shipping

Ambient Instructions

Product

Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$