

EBP11178 - IMMUNIZING PEPTIDE

SKU: EBP11178



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_000715.2; NP_963890.2; NP_963884.2; NP_963891.1; NP_963887.2; NP_963865.2; NP_963864.1; NP_963866.2; NP_001161417.1
Immunogen	Peptide with sequence ECNKQRSRHKSKDRY, from the C Terminus of the protein sequence according to NP_000715.2; NP_963890.2; NP_963884.2; NP_963891.1; NP_963887.2; NP_963865.2; NP_963864.1; NP_963866.2; NP_001161417.1.
Product Comments	The peptide was used in the production of EBP11178 - Goat Anti-CACNB2 (aa565-579) Antibody. It is currently untested in blocking applications.
Peptide Sequence	ECNKQRSRHKSKDRY
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)