

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11178 - IMMUNIZING PEPTIDE

SKU: EBP11178

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

Accession NP_000715.2; NP_963890.2; NP_963884.2; NP_963891.1; NP_963887.2; NP_963865.2; NP_963864.1;

NP_963866.2; NP_001161417.1 ID

Peptide with sequence ECNKQRSRHKSKDRY, from the C Terminus of the protein sequence according to

Immunogen NP 000715.2; NP 963890.2; NP 963884.2; NP 963891.1; NP 963887.2; NP 963865.2; NP 963864.1;

NP 963866.2; NP 001161417.1.

The peptide was used in the production of EB11178 - Goat Anti-CACNB2 (aa565-579) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

ECNKQRSRHKSKDRY Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade