

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP11207 - IMMUNIZING PEPTIDE**

**SKU:** EBP11207

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_003371.2 ID

 $\begin{tabular}{ll} \textbf{Immunogen} \\ \textbf{NP\_003371.2}. \\ \end{tabular} Peptide with sequence C-QVINETSQHHDDLE, from the C Terminus of the protein sequence according to NP\_003371.2. \\ \end{tabular}$ 

**Product** The peptide was used in the production of EB11207 - Goat Anti-Vimentin Antibody. It is currently untested in

**Comments** blocking applications.

**Peptide** 

C-QVINETSQHHDDLE Sequence

**Shipping** 

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$ Grade