

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11208 - IMMUNIZING PEPTIDE

SKU: EBP11208

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

Accession

NP 075598.2; NP 056934.2; NP 075593.1; NP 075594.1; NP 001167534.1; NP 001167535.1; NP 001167538.1 ID

Peptide with sequence RCRLRDDVQSIN, from the internal region of the protein sequence according to

Immunogen NP 075598.2; NP 056934.2; NP 075593.1; NP 075594.1; NP 001167534.1; NP 001167535.1;

NP 001167538.1.

The peptide was used in the production of EB11208 - Goat Anti-FGFR1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

RCRLRDDVQSIN Sequence

Shipping

Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade