

EBP11218 - IMMUNIZING PEPTIDE

SKU: EBP11218



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_001659.1; NP_848514.1; NP_001184254.1
Immunogen	Peptide with sequence C-SEIYHNINADQSK, from the internal region of the protein sequence according to NP_001659.1; NP_848514.1; NP_001184254.1.
Product Comments	The peptide was used in the production of EBP11218 - Goat Anti-ARNT (aa558-570) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-SEIYHNINADQSK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)