

EBP11219 - IMMUNIZING PEPTIDE

SKU: EBP11219



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001659.1;
Immunogen	Peptide with sequence CNDKERFARSDDEQS, from the internal region (near N Terminus) of the protein sequence according to NP_001659.1;.
Product Comments	The peptide was used in the production of EBP11219 - Goat Anti-ARNT (aa69-82) Antibody. It is currently untested in blocking applications.
Peptide Sequence	CNDKERFARSDDEQS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)