

EBP11250 - IMMUNIZING PEPTIDE

SKU: EBP11250



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_055148.3; NP_001091987.1
Immunogen	Peptide with sequence C-QLYTDPPQESYTT, from the internal region of the protein sequence according to NP_055148.3; NP_001091987.1.
Product	The peptide was used in the production of EBP11250 - Goat Anti-CADM1/TSLOC1 Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-QLYTDPPQESYTT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)