

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP11263 - IMMUNIZING PEPTIDE**

**SKU:** EBP11263

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_038490.2 ID

Peptide with sequence C-EARQSPPLNRFD, from the internal region (near C Terminus) of the protein sequence

according to NP\_038490.2.

The peptide was used in the production of EB11263 - Goat Anti-Adrb3 (mouse, intracellular) Antibody. It is Product

**Comments** currently untested in blocking applications.

Peptide

C-EARQSPPLNRFD Sequence

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade