



EBP11282 - IMMUNIZING PEPTIDE

SKU: EBP11282



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_006662.3
Immunogen	Peptide with sequence SLNHCTIQISELETN, from the C Terminus of the protein sequence according to NP_006662.3.
Product Comments	The peptide was used in the production of EBP11282 - Goat Anti-SLC1A7 / EAAT5 Antibody. It is currently untested in blocking applications.
Peptide Sequence	SLNHCTIQISELETN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)