

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP11306 - IMMUNIZING PEPTIDE**

**SKU:** EBP11306

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_005682.2; NP\_064693.2; NP\_001364203.1 ID

 $\label{eq:local_$ 

The peptide was used in the production of EB11306 - Goat Anti-ABCC9 / SUR2 Antibody. It is currently untested in Product

**Comments** blocking applications.

**Peptide** 

C-TSEYSINNTGKAD Sequence

Ambient Instructions Shipping

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade