

EBP11324 - IMMUNIZING PEPTIDE

SKU: EBP11324



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_002078.1

Immunogen Peptide with sequence PQDTVCDLIQSK, from the internal of the protein sequence according to NP_002078.1.

Product The peptide was used in the production of EB11324 - Goat Anti-Progranulin aa248-259 / Granulin-3 Antibody. It

Comments is currently untested in blocking applications.

Peptide Sequence PQDTVCDLIQSK

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)