

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11324 - IMMUNIZING PEPTIDE

SKU: EBP11324

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 002078.1

Immunogen Peptide with sequence PQDTVCDLIQSK, from the internal of the protein sequence according to NP_002078.1.

The peptide was used in the production of EB11324 - Goat Anti-Progranulin aa248-259 / Granulin-3 Antibody. It **Product**

Comments is currently untested in blocking applications.

Peptide

PQDTVCDLIQSK

Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png