

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP11331 - IMMUNIZING PEPTIDE**

**SKU:** EBP11331

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_002451.2; NP\_001182215.1 ID

Peptide with sequence C-DPQRQRKLITAHVEP, from the internal region (near C Terminus) of the protein sequence

according to NP\_002451.2; NP\_001182215.1.

The peptide was used in the production of EB11331 - Goat Anti-IRF4 (aa393-407) Antibody. It is currently untested **Product** 

Comments in blocking applications.

**Peptide** 

C-DPQRQRKLITAHVEP Sequence

Shipping

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$ Grade