

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11371 - IMMUNIZING PEPTIDE

SKU: EBP11371

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 036274.3

ID

Peptide with sequence C-TKASHVTPQPED, from the internal region (near C Terminus) of the protein sequence **Immunogen**

according to NP 036274.3.

The peptide was used in the production of EB11371 - Goat Anti-Cyclin D-type binding-protein 1 (aa322-333) **Product**

Comments Antibody. It is currently untested in blocking applications.

Peptide

C-TKASHVTPQPED Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade