

EBP11380 - IMMUNIZING PEPTIDE

SKU: EBP11380



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_005402.3; NP_001087239.2; NP_001158117.1; NP_001158118.1; NP_001092138.1
Immunogen	Peptide with sequence C-HLDEELQATLHDFR, from the internal region of the protein sequence according to NP_005402.3; NP_001087239.2; NP_001158117.1; NP_001158118.1; NP_001092138.1.
Product Comments	The peptide was used in the production of EBP11380 - Goat Anti-SFTPA1 / SFTPA2 (aa102-115) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-HLDEELQATLHDFR
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)