

EBP11391 - IMMUNIZING PEPTIDE

SKU: EBP11391



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_620708.1; NP_620709.1
Immunogen	Peptide with sequence C-EMVLHKVLPGRDY, from the internal region of the protein sequence according to NP_620708.1; NP_620709.1.
Product Comments	The peptide was used in the production of EB11391 - Goat Anti-MAPK9 / JNK2 beta (aa217-230) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-EMVLHKVLPGRDY
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)