

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11419 - IMMUNIZING PEPTIDE

SKU: EBP11419

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_998532.1 ID

Immunogen Peptide with sequence C-EPRQAQKHPLSAN, from the internal region of the protein sequence according to NP_998532.1.

The peptide was used in the production of EB11419 - Goat Anti-Eltd1 (Danio rerio aa298-310) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-EPRQAQKHPLSAN Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade