

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11435 - IMMUNIZING PEPTIDE

SKU: EBP11435

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_001156415.1 ID

Immunogen Peptide with sequence C-RVTVEDENDHTPT, from the internal of the protein sequence according to NP_001156415.1.

The peptide was used in the production of EB11435 - Goat Anti-Dachsous 1 / cadherin 16 (mouse) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-RVTVEDENDHTPT Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade