

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11442 - IMMUNIZING PEPTIDE

SKU: EBP11442

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_001038946.1; NP_032354.2 ID

 $\label{eq:local_$

The peptide was used in the production of EB11442 - Goat Anti-interferon-activable protein 203 (mouse isoform 1) Product

Comments Antibody. It is currently untested in blocking applications.

Peptide C-KTREQDILRHKQRP

Sequence

Ambient Instructions Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade