

# EBP11487 - IMMUNIZING PEPTIDE

SKU: EBP11487



---

## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Shipped at ambient temperature, store at -20°C.
<b>Accession ID</b>	NP_035771.1; NP_001120736.1; NP_001120734.1; NP_001120735.1; NP_001120737.1
<b>Immunogen</b>	Peptide with sequence ENNAQTQFSEPOYC, from the N Terminus of the protein sequence according to NP_035771.1; NP_001120736.1; NP_001120734.1; NP_001120735.1; NP_001120737.1.
<b>Product Comments</b>	The peptide was used in the production of EB11487 - Goat Anti-Trp63 (mouse) Antibody. It is currently untested in blocking applications.
<b>Peptide Sequence</b>	ENNAQTQFSEPOYC
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>

## DOCUMENTS

- [Data Sheet](#)