

EBP11490 - IMMUNIZING PEPTIDE

SKU: EBP11490



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_061836.2; NP_387505.1; NP_001007247.1
Immunogen	Peptide with sequence RQNCKGDSQPNKSIR, from the internal region of the protein sequence according to NP_061836.2; NP_387505.1; NP_001007247.1.
Product Comments	The peptide was used in the production of EBP11490 - Goat Anti-WDR9 / BRWD1 (aa1441-1455) Antibody. It is currently untested in blocking applications.
Peptide Sequence	RQNCKGDSQPNKSIR
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)