

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11491 - IMMUNIZING PEPTIDE

SKU: EBP11491

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_061836.2; NP_387505.1; NP_001007247.1 ID

Peptide with sequence C-NVKISETTGNSKFRP, from the internal region (near C Terminus) of the protein sequence

according to NP_061836.2; NP_387505.1; NP_001007247.1.

The peptide was used in the production of EB11491 - Goat Anti-WDR9 / BRWD1 (aa2137-2151) Antibody. It is Product

Comments currently untested in blocking applications.

Peptide

C-NVKISETTGNSKFRP Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade