

EBP11495 - IMMUNIZING PEPTIDE

SKU: EBP11495



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_033317.2
Immunogen	Peptide with sequence QVWQKAECQDPK, from the internal region of the protein sequence according to NP_033317.2.
Product Comments	The peptide was used in the production of EBP11495 - Goat Anti-Stra6 (mouse, aa187-199) Antibody. It is currently untested in blocking applications.
Peptide Sequence	QVWQKAECQDPK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)