

EBP11502 - IMMUNIZING PEPTIDE

SKU: EBP11502

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.**Unit Size** 100 µg**Storage Instructions** Shipped at ambient temperature, store at -20°C.**Accession ID** NP_116139.1**Immunogen** Peptide with sequence HPCYLDKPEEWHT, from the internal region of the protein sequence according to NP_116139.1.**Product Comments** The peptide was used in the production of EBP11502 - Goat Anti-ABHD14B (aa188-200) Antibody. It is currently untested in blocking applications.**Peptide Sequence** HPCYLDKPEEWHT**Shipping Instructions** Ambient**Product Grade** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)