

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11504 - IMMUNIZING PEPTIDE

SKU: EBP11504

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_061885.2; NP_001137293.2 ID

 $\label{eq:lower_lower_lower} \textbf{Immunogen} \ \ \frac{\text{Peptide with sequence C-RENDVKPDHETP, from the internal region of the protein sequence according to NP_061885.2; NP_001137293.2.}$

The peptide was used in the production of EB11504 - Goat Anti-PLEKHA5 / PEPP2 Antibody. It is currently untested Product

Comments in blocking applications.

Peptide

C-RENDVKPDHETP Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade