

EBP11512-B - IMMUNIZING PEPTIDE

SKU: EBP11512-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_004748.2; NP_001135888.1
Immunogen	Peptide with sequence CKQQSIAGSADSKP., from the internal region of the protein sequence according to NP_004748.2; NP_001135888.1.
Product	The peptide was used in the production of EBP11512-B - Goat Anti-AIMP1 / SCYE1 (aa137-149), Biotinylated Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	CKQQSIAGSADSKP.
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)