

EBP11514 - IMMUNIZING PEPTIDE

SKU: EBP11514



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002955.2; NP_001306125.1; NP_001306126.1; NP_001306127.1
Immunogen	Peptide with sequence C-HKYSLIKGNFH, from the internal region of the protein sequence according to NP_002955.2; NP_001306125.1; NP_001306126.1; NP_001306127.1.
Product Comments	The peptide was used in the production of EB11514 - Goat Anti-S100A8 (aa17-27) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-HKYSLIKGNFH
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png