

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP11551 - IMMUNIZING PEPTIDE**

**SKU:** EBP11551

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 081780.3; NP 001171466.1 ID

Peptide with sequence C-KSQEDTVSPTPEPQ, from the internal region (near C Terminus) of the protein **Immunogen** 

sequence according to NP 081780.3; NP 001171466.1.

The peptide was used in the production of EB11551 - Goat Anti-Prdm16 (mouse, aa1137-1150) Antibody. It is **Product** 

**Comments** currently untested in blocking applications.

**Peptide** 

C-KSQEDTVSPTPEPQ Sequence

**Shipping** 

Instructions Ambient

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png