

EBP11551 - IMMUNIZING PEPTIDE

SKU: EBP11551



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_081780.3; NP_001171466.1
Immunogen	Peptide with sequence C-KSQEDTVSPTPEPQ, from the internal region (near C Terminus) of the protein sequence according to NP_081780.3; NP_001171466.1.
Product Comments	The peptide was used in the production of EB11551 - Goat Anti-Prdm16 (mouse, aa1137-1150) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-KSQEDTVSPTPEPQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)