

# EBP11556 - IMMUNIZING PEPTIDE

**SKU:** EBP11556



---

## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage</b>	Shipped at ambient temperature, store at -20°C.
<b>Instructions</b>	
<b>Accession ID</b>	NP_005779.1; NP_001138638.1; NP_001138639.1
<b>Immunogen</b>	Peptide with sequence DPKTLISEPCGIKH, from the internal region of the protein sequence according to NP_005779.1; NP_001138638.1; NP_001138639.1.
<b>Product</b>	The peptide was used in the production of EB11556 - Goat Anti-PRMT3 (aa468-421) Antibody. It is currently untested in blocking applications.
<b>Comments</b>	
<b>Peptide Sequence</b>	DPKTLISEPCGIKH
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>

## DOCUMENTS

- [Data Sheet](#)