

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11581 - IMMUNIZING PEPTIDE

SKU: EBP11581

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 001428.1; NP 001035365.1; NP 001201831.1; NP 001201832.1 ID

Peptide with sequence C-TQDADSSRKLAH, from the internal region of the protein sequence according to Immunogen

NP 001428.1; NP 001035365.1; NP 001201831.1; NP 001201832.1.

The peptide was used in the production of EB11581 - Goat Anti-Estrogen receptor beta (aa417-428) Antibody. **Product**

Comments It is currently untested in blocking applications.

Peptide

C-TQDADSSRKLAH Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png