

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11585-B - IMMUNIZING PEPTIDE

SKU: EBP11585-B

×

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

NP_997700.1

Immunogen Peptide with sequence C-DKEFTRQPVELT., from the internal region (near C terminus) of the protein sequence

according to NP 997700.1.

Product The peptide was used in the production of EB11585-B - Goat Anti-PRKCB, Biotinylated Antibody. It is currently

Comments untested in blocking applications.

Peptide C-DKEFTRQPVELT.

Sequence Shipping

Instructions Ambient

Product

Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png