

# EBP11592 - IMMUNIZING PEPTIDE

**SKU:** EBP11592



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_065793.1

**Immunogen** Peptide with sequence C-NYGKIRGLRTPNE, from the internal region (near N Terminus) of the protein sequence according to NP\_065793.1.

**Product Comments** The peptide was used in the production of EBP11592 - Goat Anti-NLGN4X (aa51-65) Antibody. It is currently untested in blocking applications.

**Peptide Sequence** C-NYGKIRGLRTPNE

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png)

## DOCUMENTS

- [Data Sheet](#)