

EBP11605 - IMMUNIZING PEPTIDE

SKU: EBP11605



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_115674.1

Immunogen Peptide with sequence SGDYEDDLR, from the N Terminus of the protein sequence according to NP_115674.1.

Product Comments The peptide was used in the production of EBP11605 - Goat Anti-Synaptotagmin-3 Antibody. It is currently untested in blocking applications.

Peptide Sequence SGDYEDDLR

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)