

EBP11639 - IMMUNIZING PEPTIDE

SKU: EBP11639



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_001103692.1; NP_001103694.1; NP_034155.2
Immunogen	Peptide with sequence C-PKASPTPQKTS AKS, from the internal region of the protein sequence according to NP_001103692.1; NP_001103694.1; NP_034155.2.
Product Comments	The peptide was used in the production of EBP11639 - Goat Anti-Doublecortin (mouse) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-PKASPTPQKTS AKS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)