

EBP11645 - IMMUNIZING PEPTIDE

SKU: EBP11645



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001020537.2; NP_003367.4; NP_001020538.2; NP_001020539.2; NP_001028928.1; NP_001165094.1; NP_001165095.1; NP_001165096.1; NP_001165097.1; NP_001165100.1
Immunogen	Peptide with sequence KNTDSRCKARQLE, from the C Terminus of the protein sequence according to NP_001020537.2; NP_003367.4; NP_001020538.2; NP_001020539.2; NP_001028928.1; NP_001165094.1; NP_001165095.1; NP_001165096.1; NP_001165097.1; NP_001165100.1.
Product Comments	The peptide was used in the production of EBP11645 - Goat Anti-VEGFA (aa387-399) Antibody. It is currently untested in blocking applications.
Peptide Sequence	KNTDSRCKARQLE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png