

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11645 - IMMUNIZING PEPTIDE

SKU: EBP11645

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

 $NP_001020537.2; \ NP_003367.4; \ NP_001020538.2; \ NP_001020539.2; \ NP_001028928.1; \ NP_001165094.1; \ NP_001020539.2; \ NP_001028928.1; \ NP_001165094.1; \ NP_001020539.2; \ NP_001020539.$ Accession

NP_001165095.1; NP_001165096.1; NP_001165097.1; NP_001165100.1 ID

Peptide with sequence KNTDSRCKARQLE, from the C Terminus of the protein sequence according to

Immunogen NP 001020537.2; NP 003367.4; NP 001020538.2; NP 001020539.2; NP 001028928.1; NP 001165094.1;

NP 001165095.1; NP 001165096.1; NP 001165097.1; NP 001165100.1.

The peptide was used in the production of EB11645 - Goat Anti-VEGFA (aa387-399) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

KNTDSRCKARQLE Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade