

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11652 - IMMUNIZING PEPTIDE

SKU: EBP11652

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_002408.3; NP_001139438.1 ID

The peptide was used in the production of EB11652 - Goat Anti-Antigen KI-67 / MKI67 Antibody. It is currently Product

Comments untested in blocking applications.

Peptide

QRVTRSVKRCAENPK Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade