

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11653 - IMMUNIZING PEPTIDE

SKU: EBP11653

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_002408.3; NP_001139438.1 ID

 $\label{eq:local_$

The peptide was used in the production of EB11653 - Goat Anti-Antigen KI-67 / MKI67 (aa98-110) Antibody. It is Product

Comments currently untested in blocking applications.

Peptide

ENESLQNGRKSTE-C Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade